



Shown with Large Clips

RD4512



RD5512



Shown with Large Clips



RD9912

Rotators

Catalog #	Description	Voltage	RPM
099A RD4512	Motor drive & Base	120V	2-83
099A RD5512	Motor drive & Base	120V	.17-8.3
099A RD9912	Motor drive & Base	120V	2-80

Rugged Rotator

099A RD4512 & 099A RD5512

Variable speed adjusts mixing action from a gentle slosh to rapid agitation. Driven by a continuous duty motor, it can stand up to demanding schedules and long hours of operation. Some of the procedures it can be used for are tissue-culture specimens, blood cell suspensions, washing of precipitates, disintegrating tablets, extracting, dissolving slow-reacting samples and homogenizing blood solids with plasma. Containers remain visible during rotating. Sturdy, stable, well-balanced, compact base handles substantial loads. Furnished complete with the test tube head (31 clips - 10-18mm) (16 clips - 20-30mm).

- 5 lb. maximum weight capacity
- DC filtered motor for smooth operation
- Tilting motor bracket for variation of mixing action
- Variable speed control

Dry Powder Rotator - 099A RD9912

Our dry powder rotator is ideal for hard-to-mix powders. By setting the drive bar at a 20° angle, the motion of this rotator closely reproduces the mixing action achieved when mixing by hand. Not only is the powder tossed from end-to-end, it is rotated around the inside of the vessel as well. This gives a more thorough mixing action than ordinary reciprocal shakers.

- Adjustable rack fits bottles from 7" to 10" in height & up to 7" in diameter.
- 10 lb. maximum weight capacity
- DC filtered motor for smooth operation
- Silicone lined rack reduces chemical deterioration
- Tilting motor bracket for variation of mixing action
- Variable speed control

Glas-Col, LLC

711 Hulman Street P.O. Box 2128, Terre Haute, IN 47802 0128 USA

Phone: (812) 235-6167 • Fax: (812) 234-6975

Email: pinnacle@glascol.com

www.glascol.com